

Transducer Disinfection




Terason's uSmart transducers

Our family of lightweight and ergonomic transducers is ideal for easy scanning

14L3 uSmart™ Linear Array Transducer for 3200T Plus



15L4 Linear Array



15L4A Linear Array



15WL4 Linear Array



16L5 Linear Array



8L2 Linear array



4V2A Phased Array



8V3 Phased Array



5C2A Curved Array



9MC3 Curved Array



8EC4 Endocavity



16HL7 High Frequency Linear Array**



Pedoff



8TE3 Transesophageal



10EC4 Endocavity



8BP4 Endocavity



XY uSmart™ Broadband Multi-Frequency Bi-Plane



XY uSmart™ Bi-Plane Mini



Transducer Disinfection

The following products are approved for use with Terason uSmart transducers:

Ultrasound Gels*	Disinfectants*	High-Level Disinfectants**
<ul style="list-style-type: none"> Aquagel Aquasonic Clear Clear Image Clear wave Scan Ultraphonic 	<ul style="list-style-type: none"> Cavi Wipes Transeptic Spray Sono Wipes P.I. Spray II Protex Wipes Sani-Coth Wipes Sani-Cloth HB Wipes Sani-Cloth Plus Wipes Super Sani Cloth 	<ul style="list-style-type: none"> Cidex Cidex OPA Cidex Plus Metricide 14 Metricide 28 Wavicide <p>**For use only on endocavity and TEE transducers</p>

*Approved for all Terason transducers

- Trophon EPR System** has been approved for the following: 15L4, 15L4A, 16L5, 5C2A, 8V3, 9MC3, 8EC4A, 15L4, 4V2A, 8V3A, 8L2
- Trophon2** system has been approved for all probes listed for the **Trophon EPR** system, and adds 16HL7***
- When using an **enzymatic cleaner**, be sure to use the proper concentration and rinse thoroughly. Enzymatic cleaners are generically compatible for use.*
- The **Civco Astra** cleaning system can be used for endocavity probes, when used with an approved high-level disinfectant.**

***Please also see cleaning manual for 16HL7 Hockey Stick Probe Handling, Cleaning, and Sterilization for Intra-operative Use



Terason Division • Teratech Corporation
77 Terrace Hall Avenue, Burlington, MA 01803 • 1-866-TERASON • 781-270-4143 • terason.com