ACCUCARE™

REFERENCE GUIDE

EQUIPMENT:
CLASSIC STEPPER

CE

Symbol	Title of Symbol	Description of Symbol
	Manufacturer (ISO 15223-1, 5.1.1)	Indicates the medical device manufacturer.
	Authorized Representative in the European Community (ISO 15223-1, 5.1.2)	Indicates the Authorized Representative in the European Community.
	Date of manufacture (ISO 15223-1, 5.1.3)	Indicates the date when the medical device was manufactured.
	Batch code (ISO 15223-1, 5.1.5)	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog number (ISO 15223-1, 5.1.6)	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	Serial number (ISO 15223-1, 5.1.7)	Indicates the manufacturer's serial number so that a specific medical device can be identified.
	Consult instructions for use (ISO 15223-1, 5.4.3)	Indicates the need for the user to consult the instructions for use.
	European Conformity (EU MDR 2017/745, Article 20)	Indicates manufacturer declaration that the product complies with the essential requirements of the relevant European health, safety and environment protection legislation.
	Medical Device (MedTech Europe Guidance: Use of Symbols to Indicate Compliance with the MDR)	Indicates the product is a medical device.
	Quantity (IEC 60878, 2794)	To indicate the number of pieces in the package.

GENERAL EQUIPMENT INFORMATION

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

WARNING

- Before use, you should be trained in ultrasonography. For instructions on the use of your transducer, see your system's user guide.
- Prior to use, read and understand all instructions and warnings.
- Prior to use inspect device for signs of damage, if damage is evident do not use.
- Equipment is not to be modified without CIVCO authorization.
- Stepper is designed and validated for use with CIVCO accessories. For stabilizers, template grids, drapes, and other accessories, visit www.CIVCO.com.
- Stepper is packaged non-sterile and is reusable. To avoid possible patient contamination, ensure stepper is properly cleaned and disinfected before each use. See reprocessing section for instructions on how to properly clean and disinfect.
- For illustration purposes only, equipment may be shown without a drape. Always place a drape over the equipment to protect patients and users from cross-contamination.
- If the product malfunctions during use or is no longer able to achieve its intended use, stop using the product and contact CIVCO.
- Report serious incidents related to the product to CIVCO and the competent authority in your Member State, or appropriate regulatory authorities.

INTENDED USE

The equipment is intended to hold and manipulate ultrasound imaging probes, and report position, during prostate brachytherapy, cryotherapy, transperineal template-guided biopsy, and/or fiducial marker placement procedures (including volume determination of the prostate gland), and/or the application of radionuclide source(s) into the body during brachytherapy.

INDICATIONS FOR USE

Prostate - Diagnostic imaging and minimally invasive puncture procedures.

Surgical (Prostate) - Diagnostic imaging and puncture procedures.

PATIENT POPULATION

Equipment is for use in adult males with suspected or diagnosed prostate cancer.


INTENDED USERS

Equipment should be used by clinicians medically trained in ultrasound imaging. User groups may include, but are not limited to: Physicists, Radiation Oncologists, Surgeons and Urologists.

PERFORMANCE CHARACTERISTICS


- Stepper allows for incremental longitudinal and rotational movement of the ultrasound transducer to achieve its intended purpose.
- Stepper includes a quick mounting, secure connection to CIVCO stabilizers to ensure stability of the ultrasound transducer.

NOTE: For a summary of clinical benefits for this product, visit www.CIVCO.com.


- I. Grid platform
- II. Cradle latch
- III. Cradle
- IV. Cradle handle
(not included on all models)
- V. Grid rail locking knobs
- VI. Grid rail handle
- VII. Carriage scale
- VIII. Carriage
- IX. Carriage scale marker
- X. Carriage knob

Esaote TRT33
 GE Healthcare ERB, E7C8L
 Hitachi Aloka Medical C41L47RP, EUP-U533, UST-672-5/7.5
 Siemens Endo-Pil

- I. Grid platform
- II. Probe positioning locator
- III. Probe positioning notch
- IV. Cradle latch
- V. Cradle
- VI. Cradle handle
- VII. Grid rail locking knobs
- VIII. Grid rail handle
- IX. Carriage scale
- X. Carriage scale marker
- XI. Carriage knob
- XII. Carriage

BK Medical 8558, 8658, 8848, E14CL4b


- I. Grid platform
- II. Probe positioning locator
- III. Cradle latch
- IV. Cradle
- V. Grid rail locking knobs
- VI. Grid rail handle
- VII. Carriage scale
- VIII. Carriage scale marker
- IX. Carriage knob
- X. Carriage

Best® Sonalis™ TRT TriView

SETTING UP STEPPER

1. Slide stepper interface plate into stabilizer interface plate. Secure by tightening stepper quick connect knob.

- I. Stepper interface plate
- II. Stepper quick connect knob
- III. Stabilizer interface plate

2. Rotate cradle out of carriage.
3. Open the cradle latch.
4. **For Esaote TRT33, GE Healthcare ERB, E7C8L transducers:**
Place and orient transducer into the cradle, sliding transducer forward until it stops.
- For Hitachi Aloka Medical UST-672-5/7.5, Best® Sonalis™® TRT TriView transducers:**
Place and orient transducer into the cradle.
- For Hitachi Aloka Medical C41L47RP, EUP-U533, Siemens Endo-Pil transducers:**
Slide transducer back into the cradle until it stops.
- For BK Medical 8558, 8658, 8848, E14CL4b transducers:**
Slide cable end of transducer back into the cradle, ensuring locating feature fits into notch.
5. Close cradle latch and tighten knob.
6. Replace cradle with transducer back into stepper by rotating into carriage. Center transducer in rotation by feeling the center line detent during rotation.
7. Position carriage in a central location on carriage scale by turning carriage knob until carriage scale marker is at 60.
8. Loosen grid rail locking knobs and pull back on grid rail handle so grid platform is out of the way for initial transducer positioning.

DEVICE INSPECTION

1. Stepper should hold transducer securely, allow smooth rotation and precise movement along longitudinal axis of the transducer.

FUNCTIONAL PREUSE CHECKS

WARNING

- Prior to each use perform the following checks to ensure optimal stepper performance.

1. Rotate carriage knob to ensure carriage moves freely through full range of linear travel and performs stepping function.
2. Rotate cradle to ensure movement through full range of rotational travel.
3. Rotate vertical adjustment knob, grid rail locking knobs, and horizontal adjustment knobs to ensure movement through full range of rotational travel.
4. If carriage, cradle or knobs do not move freely, apply lubricant that complies with hospital policies and procedures. Super Lube® Multi-Purpose Synthetic Grease with Syncolon® (PTFE) is recommended.

- NOTE:**
- If horizontal adjustment knob does not move freely, lubricate surface areas of adjustment slide block.
 - If carriage does not move freely, lubricate outer edge of index drum.

ADJUSTING THE STEPPER AND STABILIZER

1. Grasp handles of FAM and loosen locking control knob on stabilizer.

For Micro-Touch® Stabilizer:

Adjust weight by turning weight adjustment knob on stabilizer to achieve desired feel to equipment.


2. Obtain satisfactory transducer position and starting image, tighten locking control knob on stabilizer. After stabilizer is locked into position, locking control knob should not be used until removal from the patient.
3. Use FAM to optimize position of probe.
4. Install sterile drape and grid, position grid platform appropriately and secure with grid rail locking knobs. Refer to Sterile Drape and Template Grid instructions for use.

WARNING

- Before using assembled equipment, verify needle position is correctly aligned with grid display on the ultrasound monitor. If needle position is not correctly aligned with display, adjust grid platform.
- If locking control knob on stabilizer will not hold unit securely in place, please contact CIVCO.
- Longitudinal and rotational scales are for reference only.

ADJUSTING AND SECURING GRID PLATFORM

1. Loosen Vertical Adjustment Screws with the provided hex wrench (3/32"). Align needle with the grid display on ultrasound monitor by turning the Vertical Adjustment Knob until the vertical needle position is correctly aligned with display. Tighten Vertical Adjustment Screws.
2. Loosen Horizontal Adjustment Screw with the provided hex wrench (3/32"). Align needle with the grid display on ultrasound monitor by turning the Horizontal Adjustment Knob until the horizontal needle position is correctly aligned with display. Tighten Horizontal Adjustment Screw.

- I. Vertical Adjustment Screws
- II. Vertical Adjustment Knob
- III. Horizontal Adjustment Screw
- IV. Horizontal Adjustment Knobs
- V. Adjustment Slide Block

WARNING

- Stepper is required to be cleaned and disinfected any time grid platform is readjusted and secured to ensure exposed areas of the grid platform are properly decontaminated prior to use. See Reprocessing section for instructions.

REPROCESSING

WARNING

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.
- See your system's user guide for reprocessing transducer between uses.
- These procedures for reprocessing have been validated for effectiveness and compatibility. Equipment could be damaged or cross-contaminated due to improper reprocessing.
- Do not place in a mechanical washer.
- Do not use ethylene oxide or autoclave to sterilize.

CLEANING

1. Prepare an enzymatic detergent solution, such as Enzol® Enzymatic Detergent, at manufacturer's recommendation.
2. Fully immerse Classic Stepper in prepared enzymatic detergent solution. Soak for 3 minutes.
3. After soak time, brush entire surface for 1 minute. Ensure all crevices and grooves are brushed.
4. Dispose of used enzymatic detergent solution and prepare a new enzymatic detergent solution, such as Enzol® Enzymatic Detergent, at manufacturer's recommendation.
5. Allow Classic Stepper to soak in new enzymatic detergent for 3 minutes.

6. Remove Classic Stepper from enzymatic detergent solution and rinse under running tap water for a maximum of 1 minute, but not less than 50 seconds.
7. Dry Classic Stepper with a soft, clean, lint-free cloth.

NOTE:

- If there is residual water inside the stepper, compressed air may be used at a maximum of 150 psi to remove remaining water.
- Prolonged exposure of anodized aluminum and nylon to Enzol[®] Enzymatic Detergent may cause slight discoloration of surfaces.

DISINFECTION

1. Unfold a clean germicidal wipe with at least 55% isopropyl alcohol and 0.5% quaternary ammonium, such as Super Sani-Cloth[®] Germicidal wipe, and thoroughly wet surface.
2. Allow treated surface to remain wet for no more than 2 minutes. Use additional wipe or wipes, if necessary, to ensure continuous 2 minutes wet contact time.
3. Once 2 minute wet contact time has been reached, allow surfaces to air dry.

STORAGE AND DISPOSAL

WARNING

- *Dispose of single-use components as infectious waste. Reprocess reusable components after each use according to instructions for use.*

NOTE:

- Device should be stored to keep all components together and secure.
- For questions or to order additional CIVCO products, please call +1 319-248-6757 or 1-800-445-6741 or visit www.CIVCO.com.
- All product to be returned must be in its original packaging. Contact CIVCO for further instructions as needed.

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
Esaote TRT33
 GE Healthcare ERB, E7C8L
 Hitachi Aloka Medical C41L47RP, EUP-U533, UST-672-5/7.5
 Siemens Endo-PII

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COPYRIGHT © 2020 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ACCUCARE IS A TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW.
MPS GmbH
CIVCO Medical Instruments Co., Inc.